


### catheter cushion and stabilization device



## **PIVC, MIDLINE** FOR ADULT & PEDIATRIC

### Designed to:

- Eliminate MARSI (Medical Adhesive-Related Skin Injuries)
- Prevent IV Related Pressure Injuries
- ✓ Reduce Phlebitis & Occlusion

# safe comfortable easy to use versatile

**Zpad PIV** (for most PIVs)





## **PIVC, MIDLINE** FOR ADULT & PEDIATRIC

All Zpad<sup>™</sup> adult & pediatric catheter cushions are made with state-of-the-art, skin-safe hydrocolloid adhesives to prevent skin tears and pressure injuries, and to prevent complications such as phlebitis, infiltration/extravasation and occlusion.

No other securement device can boast these unique, protective features.





#### features

- Bio-Engineered Angled Retainer
- Soft (non-latex) support
- Skin-Safe Hydrocolloid Adhesive
- Hub (blue dot) Adhesive
- Low Profile

#### benefits

- Prevents Kinking and micro-motion
- Prevents pressure injuries to skin
- Alcohol soluble to prevent skin injury
- Prevents loosening at hub catheter
- Fits beneath Transparent Dressings

### ordering information

product No	description & size range	qty
NEXJR	Nex-Zpad (BD Nexiva™)	50 pcs/box
PIVCJR	Zpad-PIV (most PIVC's)	50 pcs/box
MLJR	Zpad-pediatric (Pediatric Midlines)	50 pcs/box
ML	Zpad (Adult Midlines)	50 рся/вох

### instruction video

